

**INSTRUCTIONS FOR USE****PREFILLED SYRINGE****Product description**

Sterile, disposable prefilled syringe. Available in two options, 10ml sterile water and 10ml sterile water with 10% aqueous glycerine solution.

Intended use

The prefilled syringe provides a convenient way to inflate the balloon of Foley catheters. The prefilled syringe is beneficial in ensuring the right volume is inflated and an aseptic procedure is carried out.

REF	Description	Intended Use
20350	Prefilled Syringe w/10% Glycerine Aqueous, 10ml, Sterile	Intended for use with 100% silicone Foley catheters only. The solution reduce the permeation rate of the balloon.
20351	Prefilled Syringe w/Sterile Water, 10ml	Intended for use with latex Foley catheters only.

Shelf Life

3 years

Caution

- Do not use if the package has been opened or damaged. (sterile if pack undamaged)
- For single use only, discard after use.
- Do not use the syringe for other purposes (e.g., intramuscular and intravenous injections) except for catheter's balloon inflation.

Instruction for Use

- Wash hands and don gloves.
- Open and remove the prefilled syringe from the inner pouch.
- Remove the stopper from the luer tip. (Do not touch the tip)
- Gently, attach the luer tip to the catheter's connector.
- Then, slowly push on the thumb pad to push the plunger all the way to inflate the balloon.

Recommended disposal method

- Dispose of the used syringe into a clinical waste bag or according to the clinical waste protocol, e.g., incineration.

Recommended Storage Condition

- Store at dry, clean, and well-ventilated area at temperatures between 5° to 30°C (41°-86°F).
- Do not expose the product to radiation, sunlight, UV light, and elevated temperature as it will change the properties of the product.

Product symbol

Curas Ltd.
County Gates, Ashton Road,
Bristol BS3 2JH,
England



Manufactured by:
Curas Ltd. (RO)
69A-B Raja Uda Business
Centre Road 4, 12300
Penang, Malaysia

Ref. No.: eIFU/UR/2035X/120320/02

Date: 12-03-2020