


SPECIFICATIONS
ROMED VISITOR COATS NONWOVEN
REF: VC-L/XL/XXL



Medical device conform MDR (EU) 2017/745

Ref.no	VC-L, VC-XL, VC-XXL		
Material	Nonwoven, 35 g/m ²		
	Polypropylene 100% (Non-toxic, non-allergenic)		
Button	Plastic		
Colour	White		
Appearance	See picture: with collar, 2 pockets, 4 studs, cutting edges, no opening at back		
Physical properties fabric (MD=Machine direction, CD=Cross machine direction)	Weight	35 g/m ²	
	MD Tensile Strength	70 N/5cm	
	CD Tensile Strength	48 N/5cm	
	MD Elongation	45-95%	
	CD Elongation	45-95%	
DIMENSIONS (± 2 cm)			
(cm)	L	XL	XXL
A	100	104	108
B	64	68	72
C	58	58	60
D	28	28	29
E	8/18	8/18	8/18
F	20	20	20

*Non-woven products are not measured in thickness of the material but in weight per m²



Visitor coats can be sterilised by means of autoclaving.

These products are medical grade personal protection disposable devices and therefore are also suitable for use in the food industry.

Packing: One piece in a polybag, 50 pieces in an export carton