

# LARYVOX® TAPE COMFORT

Base Plate

- ▶ Highest possible adhesive properties
- ▶ Non-invasive
- ▶ Adapts comfortably to movements
- ▶ High flexibility and long wearing period
- ▶ Optimal for physical activity
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (KOMBI)
- ▶ Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval



 Consult instructions for use  
 CE marking  
 CE-Kennzeichnung

Product made in Germany

## Order information:

**LARYVOX® Tape COMFORT round**, pack of 15

**REF 48140**

**LARYVOX® Tape COMFORT oval**, pack of 15

**REF 48240**

**LARYVOX® Tape COMFORT XL oval**, pack of 15

**REF 48340**